

For Research Use Only

KCND3 Recombinant antibody, PBS Only (Capture)

Catalog Number:83640-2-PBS



Basic Information

Catalog Number: 83640-2-PBS	GenBank Accession Number: BC113477	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 3752	CloneNo.: 240578D2
Source: Rabbit	UNIPROT ID: Q9UK17	
Isotype: IgG	Full Name: potassium voltage-gated channel, Shal-related subfamily, member 3	
Immunogen Catalog Number: AG22183	Calculated MW: 655 aa, 73 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00616-1, KCND3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83640-2-PBS. Detection antibody: 83640-4-PBS. Standard: Ag22183. Range: 0.78-100 ng/mL