

For Research Use Only

KCND3 Recombinant antibody, PBS Only (Detector)

Catalog Number:83640-4-PBS



Basic Information

Catalog Number:

83640-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22183

GenBank Accession Number:

BC113477

GeneID (NCBI):

3752

UNIPROT ID:

Q9UK17

Full Name:

potassium voltage-gated channel,
Shal-related subfamily, member 3

Calculated MW:

655 aa, 73 kDa

Purification Method:

Protein A purification

CloneNo.:

240578E10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00616-1, KCND3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83640-2-PBS. Detection antibody: 83640-4-PBS. Standard: Ag22183. Range: 0.78-100 ng/mL