

For Research Use Only

# TADA1L Recombinant antibody, PBS Only (Detector)



Catalog Number: 83641-2-PBS

## Basic Information

Catalog Number:

83641-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14118

GenBank Accession Number:

BC015401

GeneID (NCBI):

117143

UNIPROT ID:

Q96BN2

Full Name:

transcriptional adaptor 1 (HF11 homolog, yeast)-like

Calculated MW:

335 aa, 37 kDa

Purification Method:

Protein A purification

CloneNo.:

240637H7

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00624-3, TADA1L Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83641-3-PBS. Detection antibody: 83641-2-PBS. Standard: Ag14118. Range: 0.625-80 ng/mL.