

For Research Use Only

TADA1L Recombinant antibody, PBS Only (Detector)



Catalog Number: 83641-5-PBS

Basic Information

Catalog Number: 83641-5-PBS	GenBank Accession Number: BC015401	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 117143	CloneNo.: 240637G6
Source: Rabbit	UNIPROT ID: Q96BN2	
Isotype: IgG	Full Name: transcriptional adaptor 1 (HF11 homolog, yeast)-like	
Immunogen Catalog Number: AG14118	Calculated MW: 335 aa, 37 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00624-1, TADA1L Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83641-4-PBS. Detection antibody: 83641-5-PBS. Standard: Ag14118. Range: 0.625-80 ng/mL.



Cytometric bead array standard curve of MP00624-2, TADA1L Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83641-1-PBS. Detection antibody: 83641-5-PBS. Standard: Ag14118. Range: 0.625-80 ng/mL.