

For Research Use Only

WFS1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83642-3-PBS

Basic Information

Catalog Number: 83642-3-PBS	GenBank Accession Number: BC030130	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 7466	CloneNo.: 240538B2
Source: Rabbit	UNIPROT ID: O76024	
Isotype: IgG	Full Name: Wolfram syndrome 1 (wolframin)	
Immunogen Catalog Number: AG25735	Calculated MW: 890 aa, 100 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00629-1, WFS1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83642-1-PBS. Detection antibody: 83642-3-PBS. Standard: Ag25735. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00629-3, WFS1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83642-2-PBS. Detection antibody: 83642-3-PBS. Standard: Ag25735. Range: 0.313-40 ng/mL