

For Research Use Only

HEATR3 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83643-1-PBS

Basic Information

Catalog Number: 83643-1-PBS	GenBank Accession Number: BC033077	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 55027	CloneNo.: 240561A9
Source: Rabbit	UNIPROT ID: Q7Z4Q2	
Isotype: IgG	Full Name: HEAT repeat containing 3	
Immunogen Catalog Number: AG26328	Calculated MW: 680 aa, 75 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00631-1, HEATR3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83643-1-PBS. Detection antibody: 83643-2-PBS. Standard: Ag26328. Range: 1.563-100 ng/mL



Cytometric bead array standard curve of MP00631-2, HEATR3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83643-1-PBS. Detection antibody: 83643-3-PBS. Standard: Ag26328. Range: 0.391-50 ng/mL