

For Research Use Only

PIP4K2B Recombinant antibody, PBS Only (Capture)



Catalog Number: 83647-2-PBS

Basic Information

Catalog Number:

83647-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4099

GenBank Accession Number:

BC027459

GeneID (NCBI):

8396

UNIPROT ID:

P78356

Full Name:

phosphatidylinositol-5-phosphate 4-kinase, type II, beta

Calculated MW:

416 aa, 47 kDa

Purification Method:

Protein A purification

CloneNo.:

240645C12

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00615-1, PIP4K2B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83647-2-PBS. Detection antibody: 83647-1-PBS. Standard: Ag4099. Range: 0.625-80 ng/mL.