

For Research Use Only

RNASET2 Recombinant antibody, PBS Only



Catalog Number: 83650-3-PBS

Basic Information

| | | |
|--------------------------------------|---------------------------------------|--|
| Catalog Number: 83650-3-PBS | GenBank Accession Number: BC039713 | Purification Method: Protein A purification |
| Size: 1 mg/ml | GeneID (NCBI): 8635 | CloneNo.: 240581H5 |
| Source: Rabbit | UNIPROT ID: O00584 | |
| Isotype: IgG | Full Name: ribonuclease T2 | |
| Immunogen Catalog Number: AG28353 | Calculated MW: 256 aa, 29 kDa | |

Applications

Tested Applications:
ELISA, IF/ICC

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed BxPC-3 cells using RNASET2 antibody (83650-3-RR, Clone: 240581H5) at dilution of 1:250 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) (SA00013-2). This data was developed using the same antibody clone with 83650-3-PBS in a different storage buffer formulation.