

For Research Use Only

FRK Recombinant antibody, PBS Only (Capture)



Catalog Number: 83651-1-PBS

Basic Information

Catalog Number: 83651-1-PBS	GenBank Accession Number: BC012916	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 2444	CloneNo.: 240587G6
Source: Rabbit	UNIPROT ID: P42685	
Isotype: IgG	Full Name: fyn-related kinase	
Immunogen Catalog Number: AG10372	Calculated MW: 505 aa, 56 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00622-1, FRK Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83651-1-PBS. Detection antibody: 83651-2-PBS. Standard: Ag10372. Range: 1.25-40 ng/mL