

For Research Use Only

TMEM59 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83656-1-PBS

Basic Information

| | | |
|--|---|---|
| Catalog Number: 83656-1-PBS | GenBank Accession Number: BC016374 | Purification Method: Protein A purification |
| Size: 1 mg/ml | GeneID (NCBI): 9528 | CloneNo.: 240457A7 |
| Source: Rabbit | UNIPROT ID: Q9BXS4 | |
| Isotype: IgG | Full Name: transmembrane protein 59 | |
| Immunogen Catalog Number: AG1968 | Calculated MW: 36 kDa | |

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00641-1, TMEM59 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83656-2-PBS. Detection antibody: 83656-1-PBS. Standard: Ag1968. Range: 0.78-100 ng/mL.