

For Research Use Only

LPXN Recombinant antibody, PBS Only (Detector)



Catalog Number: 83657-1-PBS

Basic Information

Catalog Number: 83657-1-PBS	GenBank Accession Number: BC019035	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 9404	CloneNo.: 240483C5
Source: Rabbit	UNIPROT ID: O60711	
Isotype: IgG	Full Name: leupaxin	
Immunogen Catalog Number: AG34956	Calculated MW: 43 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00642-1, LPXN Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83657-3-PBS. Detection antibody: 83657-1-PBS. Standard: Ag34956. Range: 0.156-20 ng/mL