

For Research Use Only

Glucagon Recombinant antibody, PBS Only (Capture)



Catalog Number: 83659-1-PBS

Basic Information

Catalog Number: 83659-1-PBS	GenBank Accession Number: BC005278	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 2641	CloneNo.: 240595B5
Source: Rabbit	UNIPROT ID: P01275	
Isotype: IgG	Full Name: glucagon	
Immunogen Catalog Number: AG8677	Calculated MW: 180 aa, 21 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00614-3, Glucagon Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83659-1-PBS. Detection antibody: 83659-4-PBS. Standard: Ag8677. Range: 0.313-40 ng/mL.