

For Research Use Only

# Glucagon Recombinant antibody, PBS Only (Detector)



Catalog Number:83659-4-PBS

## Basic Information

Catalog Number:

83659-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG8677

GenBank Accession Number:

BC005278

GeneID (NCBI):

2641

UNIPROT ID:

P01275

Full Name:

glucagon

Calculated MW:

180 aa, 21 kDa

Purification Method:

Protein A purification

CloneNo.:

240595C11

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00614-2, Glucagon Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83659-5-PBS. Detection antibody: 83659-4-PBS. Standard: Ag8677. Range: 0.313-40 ng/mL.



Cytometric bead array standard curve of MP00614-3, Glucagon Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83659-1-PBS. Detection antibody: 83659-4-PBS. Standard: Ag8677. Range: 0.313-40 ng/mL.