

For Research Use Only

MLL5 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83661-3-PBS

Basic Information

Catalog Number: 83661-3-PBS	GenBank Accession Number: BC062583	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 55904	CloneNo.: 240664A9
Source: Rabbit	UNIPROT ID: Q8IZD2	
Isotype: IgG	Full Name: myeloid/Lymphoid or mixed-Lineage leukemia 5 (trithorax homolog, Drosophila)	
Immunogen Catalog Number: AG6044	Calculated MW: 205 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00627-1, MLL5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83661-3-PBS. Detection antibody: 83661-2-PBS. Standard: Ag6044. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP00627-2, MLL5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83661-3-PBS. Detection antibody: 83661-1-PBS. Standard: Ag6044. Range: 0.625-80 ng/mL