

For Research Use Only

# FTCD Recombinant antibody, PBS Only (Capture)



Catalog Number: 83664-1-PBS

## Basic Information

Catalog Number: 83664-1-PBS	GenBank Accession Number: BC136395	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 10841	CloneNo.: 240654A6
Source: Rabbit	UNIPROT ID: O95954	
Isotype: IgG	Full Name: formiminotransferase cyclodeaminase	
Immunogen Catalog Number: AG16646	Calculated MW: 541 aa, 59 kDa	

## Applications

Tested Applications:  
Indirect ELISA, Cytometric bead array

Species Specificity:  
Human

## Background Information

### Storage

Storage:  
Store at -80°C.  
Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00636-2, FTCD Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83664-1-PBS. Detection antibody: 83664-3-PBS. Standard: Ag16646. Range: 5-80 ng/mL