

For Research Use Only

FTCD Recombinant antibody, PBS Only (Detector)



Catalog Number: 83664-3-PBS

Basic Information

Catalog Number: 83664-3-PBS	GenBank Accession Number: BC136395	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 10841	CloneNo.: 240654A11
Source: Rabbit	UNIPROT ID: O95954	
Isotype: IgG	Full Name: formiminotransferase cyclodeaminase	
Immunogen Catalog Number: AG16646	Calculated MW: 541 aa, 59 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00636-1, FTCD Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83664-2-PBS. Detection antibody: 83664-3-PBS. Standard: Ag16646. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP00636-2, FTCD Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83664-1-PBS. Detection antibody: 83664-3-PBS. Standard: Ag16646. Range: 5-80 ng/mL