

For Research Use Only

KCNH7 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83674-2-PBS

Basic Information

Catalog Number:

83674-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4532

GenBank Accession Number:

BC035815

GeneID (NCBI):

90134

UNIPROT ID:

Q9NS40

Full Name:

potassium voltage-gated channel, subfamily H (eag-related), member 7

Calculated MW:

1196 aa, 135 kDa

Purification Method:

Protein A purification

CloneNo.:

240644E8

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00649-1, KCNH7 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83674-1-PBS. Detection antibody: 83674-2-PBS. Standard: Ag4532. Range: 0.625-80 ng/mL.