

For Research Use Only

# PBRM1 Recombinant antibody

Catalog Number: 83677-2-RR



## Basic Information

Catalog Number:

83677-2-RR

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG33623

GenBank Accession Number:

BC015323

GeneID (NCBI):

55193

UNIPROT ID:

Q86U86

Full Name:

polybromo 1

Calculated MW:

180 kDa

Observed MW:

180-195 kDa

Purification Method:

Protein A purification

CloneNo.:

240655H5

Recommended Dilutions:

WB 1:2000-1:10000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse

Positive Controls:

WB : Jurkat cells, K-562 cells, HeLa cells, U2OS cells, NIH/3T3 cells

## Background Information

PBRM1, also named hPB1 BAF180 and Polybromo-1D, acts as a negative regulator of cell proliferation. PBRM1 is the second major mutated clear cell RCC cancer gene after VHL. PBRM1 loss is associated with low grade features.

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83677-2-RR (PBRM1 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.



Biolayer interferometry (BLI) kinetic assays of 83677-2-RR against Human PBRM1 were performed. The affinity constant is 2.34 nM.