

For Research Use Only

CLEC12A Recombinant antibody, PBS Only (Detector)



Catalog Number: 83680-1-PBS

Basic Information

Catalog Number: 83680-1-PBS	GenBank Accession Number: BC126291	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 160364	CloneNo.: 240791A9
Source: Rabbit	UNIPROT ID: Q5QGZ9	
Isotype: IgG	Full Name: C-type lectin domain family 12, member A	
	Calculated MW: 31 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00662-3, CLEC12A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83680-3-PBS. Detection antibody: 83680-1-PBS. Standard: Eg0024. Range: 0.625-80 ng/mL