

For Research Use Only

# FSTL3/FLRG Recombinant antibody, PBS Only (Detector)



Catalog Number: 83681-2-PBS

## Basic Information

Catalog Number:

83681-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_005860.3

GeneID (NCBI):

10272

UNIPROT ID:

O95633-1

Full Name:

follicle-stimulating-like 3 (secreted glycoprotein)

Calculated MW:

28 kDa

Purification Method:

Protein A purification

CloneNo.:

240370G1

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00667-1, FSTL3/FLRG Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83681-1-PBS. Detection antibody: 83681-2-PBS. Standard: Eg0949. Range: 0.625-80 ng/mL.