

For Research Use Only

CDCA7L Recombinant antibody, PBS Only (Capture)

Catalog Number:83682-2-PBS



Basic Information

Catalog Number: 83682-2-PBS	GenBank Accession Number: BC009352	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 55536	CloneNo.: 240639H5
Source: Rabbit	UNIPROT ID: Q96GN5	
Isotype: IgG	Full Name: cell division cycle associated 7-like	
Immunogen Catalog Number: AG10944	Calculated MW: 454 aa, 52 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00669-2, CDCA7L Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83682-2-PBS. Detection antibody: 83682-3-PBS. Standard: Ag10944. Range: 2.5-40 ng/mL.