

For Research Use Only

# B2M Recombinant antibody, PBS Only (Capture)



Catalog Number: 83683-4-PBS

## Basic Information

<b>Catalog Number:</b> 83683-4-PBS	<b>GenBank Accession Number:</b> BC032589	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 567	<b>CloneNo.:</b> 240660E11
<b>Source:</b> Rabbit	<b>ENSEMBL Gene ID:</b> ENSG00000166710	
<b>Isotype:</b> IgG	<b>UNIPROT ID:</b> P61769	
<b>Immunogen Catalog Number:</b> AG4433	<b>Full Name:</b> beta-2-microglobulin	
	<b>Calculated MW:</b> 119 aa, 14 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00673-2, B2M Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83683-4-PBS. Detection antibody: 83683-3-PBS. Standard: Ag4433. Range: 1.25-20 ng/mL