

For Research Use Only

CHI3L1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83685-3-PBS

Basic Information

Catalog Number: 83685-3-PBS	GenBank Accession Number: BC008568	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 1116	CloneNo.: 240691E3
Source: Rabbit	UNIPROT ID: P36222	
Isotype: IgG	Full Name: chitinase 3-like 1 (cartilage glycoprotein-39)	
Immunogen Catalog Number: AG33123	Calculated MW: 383 aa, 43 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00663-2, CHI3L1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83685-1-PBS. Detection antibody: 83685-3-PBS. Standard: Ag33123. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00663-3, CHI3L1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83685-2-PBS. Detection antibody: 83685-3-PBS. Standard: Ag33123. Range: 0.313-40 ng/mL