

For Research Use Only

IL6R Recombinant antibody, PBS Only (Capture)



Catalog Number:83688-1-PBS

Basic Information

Catalog Number:

83688-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC132684

GeneID (NCBI):

3570

UNIPROT ID:

P08887

Full Name:

interleukin 6 receptor

Calculated MW:

468 aa, 52 kDa

Purification Method:

Protein A purification

CloneNo.:

240734A7

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00676-2, IL6R Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83688-1-PBS. Detection antibody: 83688-2-PBS. Standard: Eg0490. Range: 0.313-40 ng/mL