

For Research Use Only

HSP47/SERPINH1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83689-1-PBS

Basic Information

Catalog Number: 83689-1-PBS	GenBank Accession Number: NM_001235.5	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 871	CloneNo.: 240752A2
Source: Rabbit	UNIPROT ID: P50454	
Isotype: IgG	Full Name: serpin peptidase inhibitor, clade H (heat shock protein 47), member 1, (collagen binding protein 1)	
	Calculated MW: 46kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00677-1, HSP47/SERPINH1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83689-2-PBS. Detection antibody: 83689-1-PBS. Standard: Eg1201. Range: 0.391-50 ng/mL