

For Research Use Only

HSP47/SERPINH1 Recombinant antibody, PBS Only (Capture)

Catalog Number:83689-4-PBS



Basic Information

Catalog Number:

83689-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001235.5

GeneID (NCBI):

871

UNIPROT ID:

P50454

Full Name:

serpin peptidase inhibitor, clade H (heat shock protein 47), member 1, (collagen binding protein 1)

Calculated MW:

46kDa

Purification Method:

Protein A purification

CloneNo.:

240752F2

Applications

Tested Applications:

WB, IHC, IF/ICC, Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human, mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Various lysates were subjected to SDS PAGE followed by western blot with 83689-4-RR (HSP47/SERPINH1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83689-4-PBS in a different storage buffer formulation.




Immunohistochemical analysis of paraffin-embedded human ovarian cancer slide using 83689-4-RR (HSP47/SERPINH1 antibody) at dilution of 1:200 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 83689-4-PBS in a different storage buffer formulation.




Immunofluorescent analysis of (4% PFA) fixed HeLa cells using HSP47/SERPINH1 antibody (83689-4-RR, Clone: 240752F2) at dilution of 1:250 and CoraLite@488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red). This data was developed using the same antibody clone with 83689-4-PBS in a different storage buffer formulation.




Immunofluorescent analysis of (4% PFA) fixed HeLa cells using HSP47/SERPINH1 antibody (83689-4-RR, Clone: 240752F2) at dilution of 1:250 and CoraLite@488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) (SA00013-2), CL594-Phalloidin (red). This data was developed using the same antibody clone with 83689-4-PBS in a different storage buffer formulation.



Sandwich ELISA standard curve of MP00677-2, Human HSP47/SERPINH1 Recombinant Matched Antibody Pair - PBS only. 83689-4-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg1201. 83689-3-PBS was HRP conjugated as the detection antibody. Range: 0.23-15 ng/mL



HepG2, HeLa and MCF-7 cells were cultured in DMEM supplemented with 10% fetal bovine serum, 100 U/mL penicillin, and 100 µg/mL streptomycin sulfate. Aliquots of the cell culture supernates were removed. Cells were lysed and assayed for human HSP47/SERPINH1. Human HSP47/SERPINH1 in HepG2 cell lysate measured 17.9 ng/mL. Human HSP47/SERPINH1 in HeLa cell lysate measured 56.9 ng/mL. Human HSP47/SERPINH1 in MCF-7 cell lysate



Biolayer interferometry (BLI) kinetic assays of 83689-4-RR against Human HSP47/SERPINH1 were performed. The affinity constant is 5.99 nM.