

For Research Use Only

# MARS Recombinant antibody, PBS Only (Capture)



Catalog Number:83690-1-PBS

## Basic Information

Catalog Number: 83690-1-PBS	GenBank Accession Number: BC002384	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 4141	CloneNo.: 240648C7
Source: Rabbit	UNIPROT ID: P56192	
Isotype: IgG	Full Name: methionyl-tRNA synthetase	
Immunogen Catalog Number: AG6619	Calculated MW: 101 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00679-3, MARS Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83690-1-PBS. Detection antibody: 83690-3-PBS. Standard: Ag6619. Range: 0.78-100 ng/mL.