

For Research Use Only

MARS Recombinant antibody, PBS Only (Detector)



Catalog Number: 83690-2-PBS

Basic Information

Catalog Number: 83690-2-PBS	GenBank Accession Number: BC002384	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 4141	CloneNo.: 240648G2
Source: Rabbit	UNIPROT ID: P56192	
Isotype: IgG	Full Name: methionyl-tRNA synthetase	
Immunogen Catalog Number: AG6619	Calculated MW: 101 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00679-2, MARS Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83690-4-PBS. Detection antibody: 83690-2-PBS. Standard: Ag6619. Range: 0.78-100 ng/mL.