

For Research Use Only

ASGR1 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83691-2-PBS

Basic Information

Catalog Number:

83691-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_001671.4

GeneID (NCBI):

432

UNIPROT ID:

P07306-1

Full Name:

asialoglycoprotein receptor 1

Calculated MW:

33 kDa

Purification Method:

Protein A purification

CloneNo.:

240672D5

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00681-1, ASGR1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83691-2-PBS. Detection antibody: 83691-3-PBS. Standard: Eg0946. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP00681-2, ASGR1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83691-2-PBS. Detection antibody: 83691-1-PBS. Standard: Eg0946. Range: 0.625-80 ng/mL