

For Research Use Only

IL18BP Recombinant antibody, PBS Only (Detector)



Catalog Number: 83692-4-PBS

Basic Information

Catalog Number: 83692-4-PBS	GenBank Accession Number: BC044215	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 10068	CloneNo.: 240711H1
Source: Rabbit	UNIPROT ID: O95998	
Isotype: IgG	Full Name: interleukin 18 binding protein	
	Calculated MW: 194 aa, 21 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00670-1, IL18BP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83692-2-PBS. Detection antibody: 83692-4-PBS. Standard: Eg1147. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP00670-2, IL18BP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83692-3-PBS. Detection antibody: 83692-4-PBS. Standard: Eg1147 Range: 0.625-80 ng/mL