

For Research Use Only

MTMR9 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 83694-1-PBS

Basic Information

Catalog Number: 83694-1-PBS	GenBank Accession Number: BC022003	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 66036	CloneNo.: 240167E6
Source: Rabbit	UNIPROT ID: Q96QG7	
Isotype: IgG	Full Name: myotubularin related protein 9	
Immunogen Catalog Number: AG34134	Calculated MW: 63 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00672-1, MTMR9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83694-2-PBS. Detection antibody: 83694-1-PBS. Standard: Ag34134. Range: 0.78-100 ng/mL



Cytometric bead array standard curve of MP00672-3, MTMR9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83694-3-PBS. Detection antibody: 83694-1-PBS. Standard: Ag34134 Range: 0.78-100 ng/mL