

For Research Use Only

MTMR9 Recombinant antibody, PBS Only (Capture)

Catalog Number: 83694-2-PBS



Basic Information

Catalog Number:

83694-2-PBS

Concentration:

1000 μ g/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34134

GenBank Accession Number:

BC022003

GeneID (NCBI):

66036

UNIPROT ID:

Q96QG7

Full Name:

myotubularin related protein 9

Calculated MW:

63 kDa

Purification Method:

Protein A purification

CloneNo.:

240167H4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Cytometric bead array standard curve of MP00672-1, MTMR9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83694-2-PBS. Detection antibody: 83694-1-PBS. Standard: Ag34134. Range: 0.78-100 ng/mL



Cytometric bead array standard curve of MP00672-2, MTMR9 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83694-2-PBS. Detection antibody: 83694-3-PBS. Standard: Ag34134. Range: 0.78-100 ng/mL



Biolayer interferometry (BLI) kinetic assay of 83694-2-PBS against Human MTMR9 was performed. The affinity constant is 0.884 nM.