

For Research Use Only

PTPN4 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83695-2-PBS

Basic Information

Catalog Number: 83695-2-PBS	GenBank Accession Number: BC010674	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 5775	CloneNo.: 240739B11
Source: Rabbit	UNIPROT ID: P29074	
Isotype: IgG	Full Name: protein tyrosine phosphatase, non-receptor type 4 (megakaryocyte)	
Immunogen Catalog Number: AG1585	Calculated MW: 106 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00680-1, PTPN4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83695-1-PBS. Detection antibody: 83695-2-PBS. Standard: Ag1585. Range: 0.156-20 ng/mL



Cytometric bead array standard curve of MP00680-3, PTPN4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83695-3-PBS. Detection antibody: 83695-2-PBS. Standard: Ag1585. Range: 0.625-20 ng/mL