

For Research Use Only

GBP4 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83699-3-PBS

Basic Information

Catalog Number: 83699-3-PBS	GenBank Accession Number: BC070055	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 115361	CloneNo.: 240730A12
Source: Rabbit	UNIPROT ID: Q96PP9	
Isotype: IgG	Full Name: guanylate binding protein 4	
Immunogen Catalog Number: AG12006	Calculated MW: 640 aa, 73 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00655-1, GBP4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83699-1-PBS. Detection antibody: 83699-3-PBS. Standard: Ag12006. Range: 0.78-50 ng/mL



Cytometric bead array standard curve of MP00655-2, GBP4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83699-2-PBS. Detection antibody: 83699-3-PBS. Standard: Ag12006. Range: 3.13-100 ng/mL