

For Research Use Only

# GBP4 Recombinant antibody, PBS Only



Catalog Number: 83699-4-PBS

## Basic Information

Catalog Number: 83699-4-PBS	GenBank Accession Number: BC070055	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 115361	CloneNo.: 240730H5
Source: Rabbit	UNIPROT ID: Q96PP9	
Isotype: IgG	Full Name: guanylate binding protein 4	
Immunogen Catalog Number: AG12006	Calculated MW: 640 aa, 73 kDa	

## Applications

Tested Applications:  
WB, ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



HeLa cells were subjected to SDS PAGE followed by western blot with 83699-4-RR (GBP4 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83699-4-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 83699-4-RR against Human GBP4 were performed. The affinity constant is 11.6 pM.