

For Research Use Only

Beta galactosidase Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: **83706-2-PBS**

Basic Information

Catalog Number: 83706-2-PBS	GenBank Accession Number: BC007493	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 2720	CloneNo.: 240788F12
Source: Rabbit	UNIPROT ID: P16278	
Isotype: IgG	Full Name: galactosidase, beta 1	
Immunogen Catalog Number: AG7792	Calculated MW: 76 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00647-1, Beta galactosidase Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83706-2-PBS. Detection antibody: 83706-3-PBS. Standard: Ag7792. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP00647-2, Beta galactosidase Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83706-5-PBS. Detection antibody: 83706-2-PBS. Standard: Ag7792. Range: 0.313-40 ng/mL