

For Research Use Only

# Mouse Crp Recombinant antibody, PBS Only (Detector)

Catalog Number:83711-2-PBS



## Basic Information

Catalog Number:

83711-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_007768

GeneID (NCBI):

12944

UNIPROT ID:

P14847

Full Name:

C-reactive protein, pentraxin-related

Purification Method:

Protein A purification

CloneNo.:

240369G12

## Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

mouse

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Sandwich ELISA standard curve of MP00815-1, mouse Crp Recombinant Matched Antibody Pair - PBS only. 83711-1-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg0756. 83711-2-PBS was HRP conjugated as the detection antibody. Range: 19.53-1250 pg/mL