

For Research Use Only

ACHE Recombinant antibody, PBS Only (Capture)



Catalog Number: 83714-1-PBS

Basic Information

Catalog Number: 83714-1-PBS	GenBank Accession Number: BC094752	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 43	CloneNo.: 240640B2
Source: Rabbit	UNIPROT ID: P22303	
Isotype: IgG	Full Name: acetylcholinesterase (Yt blood group)	
Immunogen Catalog Number: AG12146	Calculated MW: 614aa, 68 kDa; 640aa, 71 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00685-1, AChE Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83714-1-PBS. Detection antibody: 83714-2-PBS. Standard: Ag12146. Range: 1.25-80 ng/mL