

For Research Use Only

CSPG4 Recombinant antibody, PBS Only (Detector)



Catalog Number:83716-3-PBS

Basic Information

Catalog Number: 83716-3-PBS	GenBank Accession Number: NM_001897.5	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 1464	CloneNo.: 240733F2
Source: Rabbit	UNIPROT ID: Q6UVK1	
Isotype: IgG	Full Name: chondroitin sulfate proteoglycan 4	
	Calculated MW: 251kd	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00694-2, CSPG4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83716-1-PBS. Detection antibody: 83716-3-PBS. Standard: Eg0201. Range: 2.5-80 ng/mL.