

For Research Use Only

# MIB1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 83718-1-PBS

## Basic Information

Catalog Number:

83718-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2507

GenBank Accession Number:

BC022403

GeneID (NCBI):

57534

UNIPROT ID:

Q86YT6

Full Name:

mindbomb homolog 1 (Drosophila)

Calculated MW:

110 kDa

Purification Method:

Protein A purification

CloneNo.:

240676B4

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00718-3, MIB1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83718-4-PBS. Detection antibody: 83718-1-PBS. Standard: Ag2507. Range: 0.625-80 ng/mL