

For Research Use Only

MIB1 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83718-2-PBS

Basic Information

Catalog Number:

83718-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2507

GenBank Accession Number:

BC022403

GeneID (NCBI):

57534

UNIPROT ID:

Q86YT6

Full Name:

mindbomb homolog 1 (Drosophila)

Calculated MW:

110 kDa

Purification Method:

Protein A purification

CloneNo.:

240676E6

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00718-2, MIB1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83718-2-PBS. Detection antibody: 83718-4-PBS. Standard: Ag2507. Range: 1.25-80 ng/mL