

For Research Use Only

MIB1 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83718-3-PBS

Basic Information

Catalog Number: 83718-3-PBS	GenBank Accession Number: BC022403	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 57534	CloneNo.: 240676B10
Source: Rabbit	UNIPROT ID: Q86YT6	
Isotype: IgG	Full Name: mindbomb homolog 1 (Drosophila)	
Immunogen Catalog Number: AG2507	Calculated MW: 110 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00718-1, MIB1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83718-3-PBS. Detection antibody: 83718-4-PBS. Standard: Ag2507. Range: 0.625-80 ng/mL