

For Research Use Only

Mouse VCAM-1/CD106 Recombinant antibody, PBS Only (Detector)

Catalog Number: 83719-4-PBS



Basic Information

Catalog Number: 83719-4-PBS	GenBank Accession Number: CAA47989.1	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 22329	CloneNo.: 240597C2
Source: Rabbit	UNIPROT ID: P29533	
Isotype: IgG	Full Name: vascular cell adhesion molecule 1	
	Calculated MW: 81 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00701-1, MOUSE VCAM-1/CD106 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83719-1-PBS. Detection antibody: 83719-4-PBS. Standard: Eg0639. Range: 0.313-40 ng/mL.