

For Research Use Only

IL7RA/CD127 Recombinant antibody, PBS Only (Capture)



Catalog Number: 83723-3-PBS

Basic Information

Catalog Number:

83723-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

AAC83204.1

GeneID (NCBI):

3575

ENSEMBL Gene ID:

ENSG00000168685

UNIPROT ID:

P16871-1(I66T, V138I)

Full Name:

interleukin 7 receptor

Calculated MW:

52 kDa

Purification Method:

Protein A purification

CloneNo.:

240749C2

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00704-2, IL7RA/CD127 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83723-3-PBS. Detection antibody: 83723-1-PBS. Standard: Eg1279. Range: 0.625-80 ng/mL.



Cytometric bead array standard curve of MP00704-3, IL7RA/CD127 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83723-3-PBS. Detection antibody: 83723-2-PBS. Standard: Eg1279. Range: 0.625-80 ng/mL.