

For Research Use Only

# MICA Recombinant antibody, PBS Only (Detector)

Catalog Number: 83725-2-PBS



## Basic Information

Catalog Number:

83725-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

AAH16929.1

GeneID (NCBI):

100507436

Full Name:

MHC class I polypeptide-related  
sequence A

Calculated MW:

43 kDa

Purification Method:

Protein A purification

CloneNo.:

240735F10

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00686-3, MICA Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83725-1-PBS. Detection antibody: 83725-2-PBS. Standard: Eg1176. Range: 0.625-80 ng/mL