

For Research Use Only

MICA/MICB Recombinant antibody, PBS Only

Catalog Number:83725-5-PBS



Basic Information

Catalog Number:

83725-5-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

AAH16929.1

GeneID (NCBI):

4276

Full Name:

MHC class I polypeptide-related sequence A

Calculated MW:

43 kDa

Observed MW:

50-60 kDa

Purification Method:

Protein A purification

CloneNo.:

240735D9

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Background Information

Human MHC class I chain-related genes located within the HLA class I region of chromosome 6 encode MHC class I chain-related A and B (MICA and MICB) (PMID: 11429322). MICA and MICB are stress-inducible surface molecules that are not associated with β 2-microglobulin and do not present peptides (PMID: 9497295). They are expressed in intestinal epithelium and many epithelial tumors (PMID: 10359807). MICA and MICB are ligands for NKG2D which is an activating receptor that is expressed on most natural killer (NK) cells, CD8 α β T cells, and γ δ T cells (PMID: 10426993; 11491531). This antibody recognizes both MICA and MICB.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 83725-5-RR (MICA/MICB antibody) at dilution of 1:8000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 83725-5-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 83725-5-RR against Human MICA were performed. The affinity constant is 1.55 nM.