

For Research Use Only

IL1RAP Recombinant antibody, PBS Only (Capture)



Catalog Number: 83726-4-PBS

Basic Information

Catalog Number: 83726-4-PBS	GenBank Accession Number: NM_002182.4	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 3556	CloneNo.: 240756H2
Source: Rabbit	UNIPROT ID: Q9NPH3-1	
Isotype: IgG	Full Name: interleukin 1 receptor accessory protein	
	Calculated MW: 65 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00708-1, IL-1RAP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83726-4-PBS. Detection antibody: 83726-2-PBS. Standard: Eg0827. Range: 0.156-20 ng/mL



Cytometric bead array standard curve of MP00708-2, IL-1RAP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83726-4-PBS. Detection antibody: 83726-3-PBS. Standard: Eg0827. Range: 0.156-20 ng/mL