

For Research Use Only

AOX Recombinant antibody, PBS Only (Detector)



Catalog Number: 83731-2-PBS

Basic Information

| | | |
|--|--|---|
| Catalog Number: 83731-2-PBS | GenBank Accession Number: BC010425 | Purification Method: Protein A purification |
| Size: 1 mg/ml | GeneID (NCBI): 51 | CloneNo.: 240754E3 |
| Source: Rabbit | UNIPROT ID: Q15067 | |
| Isotype: IgG | Full Name: acyl-Coenzyme A oxidase 1, palmitoyl | |
| Immunogen Catalog Number: AG1405 | Calculated MW: 74 kDa | |

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00709-2, AOX Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83731-3-PBS. Detection antibody: 83731-2-PBS. Standard: Ag1405. Range: 1.56-100 ng/mL



Cytometric bead array standard curve of MP00709-3, AOX Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83731-1-PBS. Detection antibody: 83731-2-PBS. Standard: Ag1405. Range: 0.78-50 ng/mL