

For Research Use Only

AOX Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 83731-3-PBS

Basic Information

Catalog Number: 83731-3-PBS	GenBank Accession Number: BC010425	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 51	CloneNo.: 240754D4
Source: Rabbit	UNIPROT ID: Q15067	
Isotype: IgG	Full Name: acyl-Coenzyme A oxidase 1, palmitoyl	
Immunogen Catalog Number: AG1405	Calculated MW: 74 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00709-1, AOX Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83731-4-PBS. Detection antibody: 83731-3-PBS. Standard: Ag1405. Range: 0.78-100 ng/mL



Cytometric bead array standard curve of MP00709-2, AOX Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83731-3-PBS. Detection antibody: 83731-2-PBS. Standard: Ag1405. Range: 1.56-100 ng/mL