

For Research Use Only

# POT1 Recombinant antibody, PBS Only (Capture)



Catalog Number:83732-1-PBS

## Basic Information

<b>Catalog Number:</b> 83732-1-PBS	<b>GenBank Accession Number:</b> BC002923	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 1 mg/ml	<b>GeneID (NCBI):</b> 25913	<b>CloneNo.:</b> 240688B11
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q9NUX5	
<b>Isotype:</b> IgG	<b>Full Name:</b> POT1 protection of telomeres 1 homolog (S. pombe)	
<b>Immunogen Catalog Number:</b> AG0875	<b>Calculated MW:</b> 71 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00711-3, POT1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83732-1-PBS. Detection antibody: 83732-3-PBS. Standard: Ag0875. Range: 0.78-100 ng/mL